

Veterinary Anesthesia Machines ▶

R610 Tabletop Anesthesia Machine (customized)

Suitable for the inhalation anaesthesia of cats, dogs, pigs and other animals within 50 kg.



Note: Non-rebreathing circuit (for animals within 7kg)
Cycled breathing circuit (for animals above 7kg)

Order Information :

Cat No.	Product Description	Configuration	Remarks
R610IE	Veterinary Anesthesia Machine-Tabletop, Isoflurane, Easy Fil,0-4L (customized)	●	
R610IP	Veterinary Anesthesia Machine-Tabletop, Isoflurane, Pour Fil,0-4L (customized)	●	
R610SP	Veterinary Anesthesia Machine-Tabletop, Sevoflurane, Pour Fil,0-4L (customized)	●	
R610SE	Veterinary Anesthesia Machine-Tabletop, Sevoflurane, Pour Fil,0-4L (customized)	●	

Note: ● Recommended configuration ○ Optional